



Handheld 2-Function (AF-00100)

Specifications

| | |
|---------------------|--|
| Number of buttons | 4 |
| Number of functions | 2, High/Low function and leg spread function |
| Connector | 90° Hooked connector, 8-pin DIN |
| Cable | Spiralcable 8x AWG26 1285mm long |
| Materials | PPC, TPE, SS |
| Housing colour | Gray |
| Dimension LxWxH | 285x48x44mm |

Order information

| | |
|-------------------------|--|
| Article number | AF-00100 |
| Replaces Article number | HMX977/HMX977S/HMX977-1/ KPX00510.0 |
| Packaging | Per piece, 50pcs per box |
| M.O.Q. | 10 pieces |
| Deliver time | 3-5 working days |



Suitable for patient hoist

- + Sara 3000
- + Sara Plus
- + Tenor
- + Maxi Move
- + Ambulift
- + En meer..