



Handheld 1-Function (AF-00200)

Specifications

Number of buttons	2
Number of functions	1, High/Low function
Connector	Straight connector, 3-pin DIN
Cable	Spiral cable 8x AWG26 1285mm long

Materials	PPC, TPE, SS
Housing colour	Gray
Dimension LxWxH	285x48x44mm

Order information

Article number	AF-00200
Replaces Article number	KTX00890.0 KXX11600.0
Packaging	Per piece, 50pcs per box
M.O.Q.	10 pieces
Deliver time	3-5 working days



Suitable for patient hoist

+ 1 function patient hoists

