



## Handheld 2-Function (AF-00400)

### Specifications

Number of buttons	4
Number of functions	2, High/Low function and tilting function
Connector	Amphenol connector, 8-pin DIN
Cable	Spiralcable 8x AWG26 1285mm long
Materials	PPC, TPE, SS
Housing colour	Gray
Dimension LxWxH	285x48x44mm

### Order information

Article number	AF-00400
Replaces Article number	XEN145
Packaging	Per piece, 50pcs per box
M.O.Q.	10 pieces
Deliver time	3-5 working days



### Suitable for patient hoist

 Xena

