



Handheld 2-Function (AF-00500)

Specifications

Number of buttons	4
Number of functions	2, High/Low function and tilting function
Connector	Straight connector, 8-pin DIN
Cable	Spiralcable 8x AWG26 1285mm long
Materials	PPC, TPE, SS
Housing colour	Gray
Dimension LxWxH	285x48x44mm

Order information

Article number	AF-00500
Replaces	TGO.322
Article number	
Packaging	Per piece, 50pcs per box
M.O.Q.	10 pieces
Deliver time	3-5 working days



Suitable for patient hoist

- + Gemini Double
- + Older medibol hoists
- +
- +
- +
- +