



## Handheld 1-Function (AF-00600)

### Specifications

Number of buttons	2
Number of functions	1, High/Low function
Connector	90° Hooked connector, 8-pin DIN
Cable	Spiralcable 8x AWG26 1285mm long

Materials	PPC, TPE, SS
Housing colour	Gray
Dimension LxWxH	285x48x44mm

### Order information

Article number	AF-00600
Replaces Article number	8560008
Packaging	Per piece, 50pcs per box
M.O.Q.	10 pieces
Deliver time	3-5 working days



### Suitable for patient hoist

- + Carino shower chair <2012
- + Carino shower chair >2012
- +
- +
- +
- +