



## Handheld 2-Function (AF-00700)

### Specifications

Number of buttons	4
Number of functions	2, High/Low function and electrical brakes function
Connector	Hooked connector, 8-pin DIN
Cable	Spiralcable 8x AWG26 1285mm long
Materials	PPC, TPE, SS
Housing colour	Gray
Dimension LxWxH	285x48x44mm

### Order information

Article number	AF-00700
Replaces Article number	8560004
Packaging	Per piece, 50pcs per box
M.O.Q.	1 piece
Deliver time	2-5 working days



### Suitable for patient hoist

- + Alenti
- +
- +
- +
- +
- +