



Battery pack NiMH 24V 2.5Ah (AH-00300)

Specifications

Voltage	24 Volt
Capacity	2.5 Ampère/uur
Amount of cells	10 Battery cells
Type of battery	NiMH
Charging time	2.5 hours with charger C-AH-00502
Materials	ABS, Brass, PVC, Copper
Case colour	Gray RAL7014
Dimensions LxWxH	315x103x30mm

Order information

Article number	AH-00300
Compatible with	NEA0100-XX
Packaging	Individual packed, 140 pcs per pallet
M.O.Q.	1 piece
Delivery time	3-5 working days



Suitable for Patient Hoists

- + Carendo
- + Carino
- + Maxi twin
- + Maxi twin compact