



Charger 24V SLA (C-AH-00204)

Specifications

Input voltage	100-240 VAC 50-60Hz
Input current	700-500mA
Output voltage	24 Volt DC
Output current	1000mA
IP Class charger	IP40
Certification	UL2601-1 IEC60601-1 IEC60335-2-29
Operational temperature	0°C to +40°C
Storage temperature	-40°C to +70°C
Case colour	Grey RAL9006, Black charger
Dimensions LxWxH	715x100x45mm



Order Information

Article number	C-AH-00200 for EU C-AH-00202 for UK C-AH-00204 for USA C-AH-00203 for AUS
Compatible with	NDA8200
Packaging	Individual packed, 140 pcs per pallet
M.O.Q.	10 pieces
Delivery time	3-5 working days

Features

- + Suitable for AH-00100 / NDA0100-XX HMX991
- + Charger specially designed for medical equipment
- + Delivered with connection cable suited for UK, USA, EU or AUS
- + When the battery pack is charged, the charger will continue as a trickle charger