



Charger 24V SLA (C-AH-00404)

Specifications

Input voltage	100-240 VAC 50-60Hz
Input current	700-500mA
Output voltage	24 Volt DC
Output current	1000mA
IP Class charger	IP40
Certification	UL2601-1 IEC60601-1 IEC60335-2-29
Operational temperature	0°C to +40°C
Storage temperature	-40°C to +70°C
Case colour	White RAL9010, Black charger
Dimensions LxWxH	330x95x80mm



Order Information

Article number	C-AH-00400 for EU C-AH-00402 for UK C-AH-00404 for USA C-AH-00403 for AUS
Compatible with	KTA0101
Packaging	Individual packed, 100 pcs per pallet
M.O.Q.	10 pieces
Delivery time	3-5 working days

Features

- + Suitable for AH-00700 / KTA0102
- + Charger specially designed for medical equipment
- + Delivered with connection cable suited for UK, USA, EU or AUS
- + When the battery pack is charged, the charger will continue as a trickle charger