



## Sling Scale Patient Hoist (sc-00150)

### Specifications

Accuracy Class according to OIML R76	C2 / C1.6
Maximum capacity	C2 = 200KG / C1.6 = 320KG
Maximum load	350KG
Breaking load	1000KG
Minimum interval	C2 = 100gr / C1.6 = 200gr
Display	5 digits, digital
Battery	4x AA batteries
Battery Life	270 hours (1600mA batteries)
Battery Voltage	3.6—6V
Current consumption activated	< 6mA
Quiescent consumption standby	< 0.001mA
Protection DIN40050	IP54
Operating temp.	-20°C tot +60°C
Storage temp.	-40°C tot +85°C
Dimensions LxWxH	165x70x116mm

### Order information

Article number	SC-00150
Replaces Art. number	
Packaging	Individually packed, 10pcs per box
M.O.Q.	1 pcs
Delivery Time	2 weeks

### Accessories

Adapter M12 threaded	SC-00201
Adapter Yoke 14mm	SC-00202
Adapter hook	SC-00203
Adapter eye	SC-00204
Calibration	SC-00206



### Suited for patient hoists

**+** All hoists with scale support

